

CIAAlign - Clean and Interpret Alignments

Charlotte Tumescheit, Dr. Andrew E. Firth, Dr. Katherine Brown

CIAAlign is a command line tool which performs various functions to clean and analyse a multiple sequence alignment (MSA).

Sign up here for updates when a new feature is added to CIAAlign

The tool is designed to be highly customisable, allowing users to specify exactly which functions to run and which settings to use. It is also transparent, generating a clear log file and alignment markup showing exactly how the alignment has changed and what has been removed by which function.

This allows the user to:

- Remove sources of noise from their MSA
 - Remove insertions which are not present in the majority of sequences
 - Remove sequences below a threshold number of bases or amino acids
 - Crop poorly aligned sequence ends
 - Remove columns containing only gaps
 - Remove sequences above a threshold level percentage of divergence from the majority
 - Remove either end of an alignment where columns don't meet a minimum identity threshold and coverage level
- Generate consensus sequences
- Visualise alignments
 - Generate image files showing the alignment before and after using CIAAlign cleaning functions and showing which columns and rows have been removed
 - Draw sequence logos
 - Visualise coverage at each position in the alignment
- Analyse alignment statistics
 - Generate a similarity matrix showing the percentage identity between each sequence pair
- Unalign the alignment
- Replace U's by T's

Citation

If you found CIAAlign useful, please cite:

Tumescheit C, Firth AE, Brown K. 2022. CIAAlign: A highly customisable command line tool to clean, interpret and visualise multiple sequence alignments. PeerJ 10:e12983 <https://doi.org/10.7717/peerj.12983>

Requirements

- python \geq 3.6
- matplotlib \geq 2.1.1
- numpy \geq 1.16.3
- scipy \geq 1.3.0

Installation

The easiest way to install CIAAlign is using conda or pip3.

Conda

```
conda install -c bioconda cialign
```

[link](#)

pip3

```
pip3 install cialign
```

[link](#)

The current release of CIAAlign can also be downloaded directly using this [link](#),

If you download the package directly, you will also need to add the CIAAlign directory to your PATH environment variable as described [here](#)

Usage

Basic Usage

```
CIAAlign --infile INFILE --outfile_stem STEM --inifile my_config.ini
```

Parameters Parameters can be specified in the command line or in a config file using the naming system below.

A template config file is provided in CIAAlign/templates/ini_template.txt - edit this file and provide the path to the --infile argument. If this argument is not provided command line arguments and defaults will be used.

Parameters passed in the command line will take precedence over config file parameters, which take precedence over defaults.

Command help can be accessed by typing `CIAAlign --help`

Parameter	Description	Default
--infile	Path to input alignment file in FASTA format	None
--inifile	Path to config file	None
--outfile_stem	Prefix for output files, including the path to the output directory	CIAAlign
--silent	Do not print progress to the screen	False
--all	Use all available functions with default parameters. Does not currently include crop_divergent	False
--clean	Use all available cleaning functions (except crop_divergent) with default parameters	False
--visualise	Use all available mini alignment visualisation functions with default parameters	False
--interpret	Use all available interpretation functions (except sequence logos) with default parameters	False
--help	Show all available parameters with an explanation	None
--version	Show the version	None

Beside these main parameters, the use of every function and corresponding thresholds can be specified by the user by adding parameters to the command line or by setting them in the configuration file. Available functions and their parameters will be specified in the following section.

CIAAlign always produces a log file, specifying which functions have been run with which parameters and what has been removed. It also outputs a file that only specifies what has been removed with the original column positions and the sequence names.

Output files:

- **OUTFILE_STEM_log.txt** - general log file
- **OUTFILE_STEM_removed.txt** - removed columns positions and sequence names text file

Cleaning an MSA

Each of these steps (if specified) will be performed sequentially in the order specified in the table below.

The “cleaned” alignment after all steps have been performed will be saved as **OUTFILE_STEM_cleaned.fasta**

remove_divergent, remove_insertions, crop_ends and crop_divergent require three or more sequences in the alignment, remove_short and remove_gap_only require two or more sequences.

Parameter	Description	Default Value	Min	Max
--remove_divergent	Remove sequences with \leq remove_divergent_minperc positions at which the most common base / amino acid in the alignment is present	False	NA	NA
--remove_divergent_minperc	Minimum proportion of positions which should be identical to the most common base / amino acid in order to be preserved	0.65	0	1
--remove_insertions	Remove insertions found in \leq insertion_min_perc of sequences from the alignment	False	NA	NA
--remove_divergent_retain	Do not remove sequences with this name when running the remove divergent function	None	NA	NA
--remove_divergent_retain_str	Do not remove sequences with names containing this character string when running the remove divergent function	None	NA	NA
--remove_divergent_retain_list	Do not remove sequences with names listed in this file when running the remove divergent function	None	NA	NA
--remove_insertions	Remove insertions found in \leq insertion_min_perc of sequences from the alignment	False	NA	NA
--insertion_min_size	Only remove insertions \geq this number of residues	3	1	n_col
--insertion_max_size	Only remove insertions \leq this number of residues	200	1	1000
--insertion_min_flank	Minimum number of bases on either side of an insertion to classify it as an insertion	5	0	n_col/2
--insertion_min_perc	Remove insertions which are present in less than this proportion of sequences	0.5	0	1
--crop_ends	Crop the ends of sequences if they are poorly aligned	False	NA	NA
--crop_ends_mingap_perc	Minimum proportion of the sequence length (excluding gaps) that is the threshold for change in gap numbers.	0.05	0	0.6
--crop_ends_redefine_perc	Proportion of the sequence length (excluding gaps) that is being checked for change in gap numbers to redefine start/end.	0.1	0	0.5
--crop_ends_retain	Do not crop sequences with this name when running the crop ends function	None	NA	NA
--crop_ends_retain_str	Do not crop sequences with names containing this character string when running the crop ends function	None	NA	NA
--crop_ends_retain_list	Do not crop sequences with names listed in this file when running the crop ends function	None	NA	NA
--crop_divergent	Crop either end of the alignment until $>$ crop_divergent_min_prop_ident residues in a column are identical and $>$ crop_divergent_min_prop_nongap residues are not gaps, over buffer_size consecutive columns	False	NA	NA
--crop_divergent_min_prop_ident	Minimum proportion of identical residues in a column to be retained by crop_divergent	0.5	0.01	1
--crop_divergent_min_prop_nongap	Minimum proportion of non gap residues in a column to be retained by crop_divergent	0.5	0.01	1
--crop_divergent_buffer_size	Minimum number of consecutive columns which must meet the criteria for crop_divergent to be retained	5	1	n_col
--remove_short	Remove sequences \leq remove_min_length amino acids from the alignment	False	NA	NA
--remove_min_length	Sequences are removed if they are shorter than this minimum length, excluding gaps.	50	0	n_col
--remove_short_retain	Do not remove sequences with this name when running the remove short function	None	NA	NA
--remove_short_retain_str	Do not remove sequences with names containing this character string when running the remove short function	None	NA	NA

Parameter	Description	Default Value	Min	Max
<code>--remove_short_retain_list</code>	Do not remove sequences with names listed in this file when running the remove short function	None	NA	NA
<code>--keep_gaponly</code>	Keep gap only columns in the alignment	False	NA	NA
<code>--retain</code>	Do not edit or remove sequences with this name when running any rowwise function (currently remove divergent, crop ends and remove short)	None	NA	NA
<code>--retain_str</code>	Do not edit or remove sequences with names containing this character string when running any rowwise function	None	NA	NA
<code>--retain_list</code>	Do not edit or remove sequences with names listed in this file when running any rowwise function	None	NA	NA

Note: if the sequences are short (e.g. < 100), a low `crop_ends_mingap_perc` (e.g. 0.01) will result in a change of gap numbers that is too low (e.g. 0). If this happens, the change in gap numbers will be set to 2 and a warning will be printed.

Generating a Consensus Sequence

This step generates a consensus sequence based on the cleaned alignment. If no cleaning functions are performed, the consensus will be based on the input alignment. For the “majority” based consensus sequences, where the two most frequent characters are equally common a random character is selected.

Output files:

- `OUTFILE_STEM_consensus.fasta` - the consensus sequence only
- `OUTFILE_STEM_with_consensus.fasta` - the cleaned alignment plus the consensus

Parameter	Description	Default
<code>--make_consensus</code>	Make a consensus sequence based on the cleaned alignment	False
<code>--consensus_type</code>	Type of consensus sequence to make - can be majority, to use the most common character at each position in the consensus, even if this is a gap, or majority_nongap, to use the most common non-gap character at each position	majority
<code>--consensus_keep_gaps</code>	If there are gaps in the consensus (if majority_nongap is used as consensus_type), should these be included in the consensus (True) or should this position in the consensus be deleted (False)	False
<code>--consensus_name</code>	Name to use for the consensus sequence in the output fasta file	consensus

Unaligning the Alignment

This function simply removes the gaps from the input or output alignment and creates an unaligned file of the sequences.

Output files:

- `OUTFILE_STEM_unaligned_input.fasta` - unaligned sequences of input alignment
- `OUTFILE_STEM_unaligned_output.fasta` - unaligned sequences of output alignment

Parameter	Description	Default
<code>--unalign_input</code>	Generates a copy of the input alignment with no gaps	False
<code>--unalign_output</code>	Generates a copy of the output alignment with no gaps	False

Replacing U's by T's

This function replaces the U nucleotides by T nucleotides without disturbing the sequence names.

Output files:

- `OUTFILE_STEM_T_input.fasta` - input alignment with T's instead of U's
- `OUTFILE_STEM_T_output.fasta` - output alignment with T's instead of U's

Parameter	Description	Default
<code>--replace_input</code>	Generates a copy of the input alignment with T's instead of U's	False
<code>--replace_output</code>	Generates a copy of the output alignment with T's instead of U's	False

Visualising Alignments

Each of these functions produces some kind of visualisation of your alignment.

Mini Alignments

These functions produce “mini alignments” - images showing a small representation of your whole alignment, so that gaps and poorly aligned regions are clearly visible.

Output files:

- `OUTFILE_STEM_input.png` (or `svg`, `tiff`, `jpg`) - the input alignment
- `OUTFILE_STEM_output.png` (or `svg`, `tiff`, `jpg`) - the cleaned output alignment
- `OUTFILE_STEM_markup.png` (or `svg`, `tiff`, `jpg`) - the input alignment with deleted rows and columns marked

Parameter	Description	Default
<code>--plot_input</code>	Plot a mini alignment - an image representing the input alignment	False
<code>--plot_output</code>	Plot a mini alignment - an image representing the output alignment	False
<code>--plot_markup</code>	Draws the input alignment but with the columns and rows which have been removed by each function marked up in corresponding colours	False
<code>--plot_dpi</code>	DPI for mini alignments	300
<code>--plot_format</code>	Image format for mini alignments - can be png, svg, tiff or jpg	png
<code>--plot_width</code>	Mini alignment width in inches	5
<code>--plot_height</code>	Mini alignment height in inches	3
<code>--plot_keep_numbers</code>	Label rows in mini alignments based on input alignment, rather than renumbering	False
<code>--plot_force_numbers</code>	Force all rows in mini alignments to be numbered rather than labelling e.g. every 10th row for larger plots Will cause labels to overlap in larger plots	False

Sequence logos

These functions draw sequence logos representing your output (cleaned) alignment. If no cleaning functions are specified, the logo will be based on your input alignment.

Output_files:

- `OUTFILE_STEM_logo_bar.png` (or `svg`, `tiff`, `jpg`) - the alignment represented as a bar chart
- `OUTFILE_STEM_logo_text.png` (or `svg`, `tiff`, `jpg`) - the alignment represented as a standard sequence logo using text

Parameter	Description	Default
<code>--make_sequence_logo</code>	Draw a sequence logo	False
<code>--logo_start</code>	Start of sequence logo	0
<code>--logo_end</code>	End of sequence logo	MSA length
<code>--sequence_logo_type</code>	Type of sequence logo - bar/text/both	bar
<code>--sequence_logo_dpi</code>	DPI for sequence logo	300
<code>--sequence_logo_font</code>	Font (see NB below) for bases / amino acids in a text based sequence logo	monospace
<code>--sequence_logo_nt_per_row</code>	Number of bases / amino acids to show per row in the sequence logo, where the logo is too large to show on a single line	50
<code>--sequence_logo_filetype</code>	Image file type to use for the sequence logo - can be png, svg, tiff or jpg	png

NB: to see available fonts on your system, run `CIAAlign -list_fonts_only` and view `CIAAlign_fonts.png`

Coverage Plots

This function plots the number of non-gap residues at each position in the alignment.

Output file:

- `OUTFILE_STEM_input_coverage.png` (or `svg`, `tiff`, `jpg`) - image showing the input alignment coverage
- `OUTFILE_STEM_output_coverage.png` (or `svg`, `tiff`, `jpg`) - image showing the output alignment coverage

Parameter	Description	Default
<code>--plot_coverage_input</code>	Plot the coverage of the input MSA	False
<code>--plot_coverage_output</code>	Plot the coverage of the output MSA	False
<code>-plot_coverage_dpi</code>	DPI for coverage plot	300
<code>-plot_coverage_height</code>	Height for coverage plot (inches)	3
<code>-plot_coverage_width</code>	Width for coverage plot (inches)	5
<code>-plot_coverage_colour</code>	Colour for coverage plot (hex code or name)	#007bf5
<code>-plot_coverage_filetype</code>	File type for coverage plot (png, svg, tiff, jpg)	png

Analysing Alignment Statistics

These functions provide additional analyses you may wish to perform on your alignment.

Similarity Matrices

Generates a matrix showing the proportion of identical bases / amino acids between each pair of sequences in the MSA.

Output file:

- `OUTFILE_STEM_input_similarity.tsv` - similarity matrix for the input file
- `OUTFILE_STEM_output_similarity.tsv` - similarity matrix for the output file

Parameter	Description	Default
<code>--make_similarity_matrix_input</code>	Make a similarity matrix for the input alignment	False
<code>--make_similarity_matrix_output</code>	Make a similarity matrix for the output alignment	False
<code>-make_simmatrix_keepgaps</code>	0 - exclude positions which are gaps in either or both sequences from similarity calculations, 1 - exclude positions which are gaps in both sequences, 2 - include all positions	0
<code>-make_simmatrix_dp</code>	Number of decimal places to display in the similarity matrix output file	4
<code>-make_simmatrix_minoverlap</code>	Minimum overlap between two sequences to have non-zero similarity in the similarity matrix	1